



# IVF MEDIA PRODUCTS

*... for the Next Generation*

## InVitroCare® IVF MEDIA SERIES

### HTF HEPES

For **egg retrieval, follicle flushing, ICSI procedures and cryopreservation**

Formulated with EDTA and a stabilized form of glutamine

12 month shelf life at 2-8°C

Cat # 2002 - 100ml or Cat # 2002-5 - 500ml

### HTF

For use as a **fertilization medium** and for early embryo culture through 2PN stage

Formulated with EDTA and a stabilized form of glutamine

120 day shelf life at 2-8°C

Cat # 2001 - 100ml

### IVC-ONE™

For **early embryo development** from Day 0 through compaction

**Glucose and Phosphate Free**

Formulated with EDTA, Glutathione, and a stabilized form of glutamine

120 day shelf life at 2-8°C

Cat # 2006 - 100ml

### IVC-TWO™

For use in **Day 3 embryo transfers** and ICSI procedures

**Low Glucose and Phosphate Free**

Formulated with EDTA, Glutathione, and a stabilized form of glutamine

120 day shelf life at 2-8°C

Cat # 2008 - 100ml

### IVC-THREE™

For **extended embryo culture** and transfer through expanded **Blastocyst (Day 5)**

**Phosphate Free and elevated levels of Glucose**

Supplemented with amino acids and vitamins

Formulated with Glutathione

120 day shelf life at 2-8°C

Cat # 2007 - 100ml

8415 Progress Drive, Suite F  
Frederick, MD 21701  
800 786 1986 ~ 301 360 9660  
fax 301 360 9661  
[www.invitrocare.com](http://www.invitrocare.com)