

IUI MADE SIMPLE...

... Laboratory Results in an Office Setting

InVitroCare[®] SpermCare[™] One-Step

Catalog # 2224

- ◆ A *ready-to-use* sperm preparation system for INTRAUTERINE INSEMINATION (IUI)
- ◆ Includes *pre-measured* gradient in 15ml centrifuge tubes and *pre-measured* InVitroCare[®] Sperm Wash Medium in glass vials

ADD SEMEN

SPIN

WASH

LOAD CATHETER



- ◆ Each SpermCare[™] One-Step System contains enough product for 5 determinations.
- ◆ SpermCare[™] One-Step is designed to reduce the chance of contamination and helps eliminate product waste.
- ◆ Gradient, consisting of silane-coated colloidal silica particles suspended in HEPES-buffered HTF, removes bacteria, epithelial cells, cell debris, immature sperm and lymphocytes
- ◆ All system components contain glucose as a usable energy substrate to maintain normal cell metabolism and sperm function
- ◆ Aseptically packaged and terminally sterilized

This product, and other IUI media products, catheters, and disposables are available online at www.iuimadesimple.com

8415 Progress Drive, Suite F
Frederick, MD 21701
800 786 1986 ~ 301 360 9660
fax 301 360 9661
www.invitrocare.com